

REF 43-32

**For cleaning and disinfection  
of hard surfaces of non-invasive,  
non-critical and non-lumened  
medical devices. Effective against  
bacteria, yeast and fungi.**

#### Intended Use

Protex® Pro Disinfectant Spray is a solution for cleaning and disinfecting hard surfaces of non-invasive, non-critical and non-lumened medical devices.

#### Applications

The applications include – but are not limited to – skin surface ultrasound transducers and stations, IV poles and other hard non-porous, non-critical surfaces.

This product is not intended to be used to disinfect invasive medical devices.

#### Contents

##### Active ingredients:

Didecyl dimethyl ammonium chloride... 0.23%

Alkyldimethyl benzyl ammonium chloride... 0.13%

Polyhexamethylene biguanide... 0.208%



500ml RTU

A graphic of a medical device, possibly a transducer, with three lines radiating from its top, positioned above the 'X' in the brand name.  
**protex<sup>®</sup> pro**  
disinfectant spray

LA-043-0001A



**Directions for Use** Wear gloves when handling disinfectants and medical devices. • Remove any solid waste from the surface. • Apply the product to the surface and spread with a low-linting wipe, ensuring all areas are wet for the five-minute contact time. Leave surface to dry. • Dispose wipe to clinical waste – do not macerate, do not flush, do not reuse. • For professional use only.

**Precautions** Wear protective gloves, protective clothing and eye protection. • Do not use on live electrical components. • Dispose in accordance with local policy and national regulations. • Refer to the Safety Data Sheet of the product, as applicable.

**Biocidal Performance** Protex<sup>®</sup> Pro Disinfectant Spray is effective against the following microorganisms of concern and also destroys DNA and RNA. • *Candida albicans* • *Trichophyton interdigitale* • *Acinetobacter baumannii* • *Campylobacter jejuni* • *Enterococcus faecalis* • *Enterococcus hirae* • *Escherichia coli* • *Klebsiella pneumoniae* • Methicillin-resistant *Staphylococcus aureus* (MRSA) • Methicillin-resistant *Staphylococcus pseudintermedius* (MRSP) • *Proteus vulgaris* • *Pseudomonas aeruginosa* • *Rhodococcus equi* • *Salmonella typhimurium* • *Staphylococcus aureus* • *Stenotrophomonas maltophilia* • *Streptococcus equi* • Vancomycin-resistant *Enterococci* (VRE) *Enterococcus faecium*

#### Warning

- Contains polymeric biguanide hydrochloride.
- May produce an allergic reaction.
- Causes serious eye irritation.

#### Precautionary statements

- Wear protective gloves/protective clothing/eye protection.
- **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do-continue rinsing.
- If eye irritation persists get medical advice.



2797



**Parker Laboratories, Inc.**, 286 Eldridge Road, Fairfield, New Jersey 07004 • USA



Bedrijvenpark Twente 305, 7602 KL Almelo, The Netherlands