

**The GROUNDBREAKING
Sterile Barrier and
Securement Dressings
UNIQUELY DESIGNED for UGPIV**



ULTRA DRAPE®
UGPIV Barrier and Securement

ULTRA DRAPE® II
UGPIV Barrier and Securement



ULTRA DRAPE®

UGPIV Barrier and Securement

The Dressing That Simplifies UGPIV

UltraDrape was engineered to facilitate a user-friendly, no-touch aseptic procedure while...

- ▶ Offering a streamlined approach to barrier and dressing securement
- ▶ Reducing procedure time and cost when compared with use of a sterile gel and sterile cover
- ▶ Eliminating time-consuming clean up
- ▶ Minimizing dressing securement failure



ULTRA DRAPE® II

UGPIV Barrier and Securement

The Evolutionary Dressing for UGPIV

Reengineered to improve visibility, versatility, and functionality...

- ▶ Smaller size for use on a broader patient population
- ▶ Flex-Tab™ reinforcement tabs improve “standability” creating a secure barrier between transducer and skin
- ▶ Borderless end for ease of scanning in both long and short axis. Allows for visualization of “walking the catheter” into the vessel and imaging beyond the dressing.



EASY AS 1, 2, 3!



1

Mark the selected site and adjust the gain brighter. Peel and fold off tab #1.




2

Position fold edge of UltraDrape on mark and stick to skin. Apply gel to back area, tab #2.



3

Insert, peel gel layer off and pull tab #3 down.
FINISHED!



STRONGLY AGREED

The barrier dressing provided gel and probe separation from the skin



**87%
STRONGLY
AGREED**

UltraDrape improved patient care by facilitating aseptic technique

- ▶ **Speeds UGPIV insertion with full safety in mind**
- ▶ **Promotes achievable standardization of UGPIV insertions for all departments**
- ▶ **Sterile protection for transducer enabling low-level disinfection**
- ▶ **Safe and cost-effective for all**
- ▶ **Easily added to vendor-prepared, customized IV kits**



PROCEDURAL VIDEO

parkerlabs.com/products/ultradrape



SUPPORT/TRAINING

To request additional samples, a demo or training support contact


customerservice@parkerlabs.com



COST-SAVINGS CALCULATOR

parkerlabs.com/products/ultradrape

Do you prefer using a sterile transducer cover or UltraDrape?



Reference: Drafz M, Goeller K, Dizon B, Cobbs A, Moureau N (2023) Multi-center Qualitative Observational Evaluation of Ultrasound Probe Protection using a Sterile Transparent Barrier and Securement Dressing to Standardize UGPIV Catheter Insertions. Int J Nurs Health Care Res 6: 1418. <https://protect-us.mimecast.com/s/Aa89CJ6kOmH8O0ZUVQKm6?do-main=doi.org> <https://doi.org/10.29011/2688-9501.101418>

The uniquely designed, dual-action barrier and securement dressings for use during UGPIV procedures. Dressings provide sterile probe protection and keep gel away from the insertion site.

ULTRA DRAPE®

UGPIV Barrier and Securement



Product Number

34-15

Package Size

5.8 in x 3.25 in (14.72 cm x 8.25 cm)
Box of 50, 4 bx/cs

ULTRA DRAPE® II

UGPIV Barrier and Securement



Product Number

34-10

Package Size

4.68 in x 2.63 in (11.9 cm x 6.7 cm)
Box of 50, 4 bx/cs



parkerlabs.com/ultradrape

parkerlabs.com/ultradrape2



Parker Laboratories, Inc.

The sound choice in patient care.™