
UltraDrape® and UltraDrape® II UGPIV Barrier and Securement

REF 34-15, 34-10

- | | |
|---------------------------------------|-----------------------------------|
| EN Instructions for use | JA 使用方法 |
| ES Instrucciones de uso | KO 적용 사항 |
| DE Gebrauchsanweisung | ZH 使用说明 |
| FR Mode d'emploi | AR إرشادات الاستعمال |
| IT Istruzioni per l'uso | BG Показания за употреба |
| PT Instruções de utilização | HR Indikacije za uporabu |
| DA Brugsanvisning | CS Indikace k použití |
| FI Käyttöohjeet | LV Lietošanas indikācijas |
| NL Gebruiksaanwijzing | LT Naudojimo indikacijos |
| SV Indikationer för användning | RO Indicații de utilizare |
| EL Οδηγίες Χρήσης | SK Indikácie na použitie |
| HU Használati útmutatás | NO Indikasjoner for bruk |
| PL Wskazania do stosowania | TR Kullanım Endikasyonları |
| HE הוראות שימוש | |

PARKER

PARKER LABORATORIES, INC.

Parker Laboratories, Inc.
286 Eldridge Road
Fairfield, NJ 07004 USA
(973) 276-9500
parkerlabs.com

CE 2797

Patent parkerlabs.com/ultradrape

Parker Laboratories BV
Bedrijvenpark Twente 305
7602 KL Almelo
The Netherlands

Instructions for Use:**Description**

UltraDrape / UltraDrape II is a uniquely designed, sterile, dual-action barrier and securement dressing for use during Ultrasound-Guided Peripheral Intravenous (UGPIV) procedures. UltraDrape is ASAP (The Association for Safe Aseptic Practice) approved for facilitating Standard - ANTT® (Aseptic Non-Touch Technique). UltraDrape is non-irritating, non-sensitizing and is not made with natural rubber latex.

UltraDrape is 5.8 in (14.72 cm) x 3.25 in (8.25 cm)

UltraDrape II is 4.68 in (11.9 cm) x 2.63 in (6.7 cm)

Ultrasound gel is applied to removable film layer to prevent contamination of injection site. Convenient adhesive strips on the release film provide added securement of tubing and catheter hub surround. The bifurcated, "stand-alone" design prevents gel from reaching the IV site while enabling seamless catheter insertion and securement.

Indications for Use

UltraDrape / UltraDrape II is a sterile, single-use disposable device, serves as a sterile, viral barrier when used in conjunction with ultrasound transducers to shield patients from the potential transfer of contaminants and biological risk agents that may be introduced during ultrasound scanning and needle-guided procedures. It also serves as a transparent dressing to protect and secure devices to the patient's skin when required. UltraDrape / UltraDrape II is indicated for adults and pediatric patients. For use by healthcare professionals appropriately trained in ultrasound within healthcare settings utilizing ultrasound.

Any serious incident that has occurred in relation to device should be reported to the manufacturer and competent authority of the EU Member State in which user and/or patient is established.

**Precautions**

- Discard product if package is damaged or evidence of contamination is present.
- Do not re-use.
- Do not re-sterilize.

UltraDrape / UltraDrape II Application Considerations

- Assure skin is clean, free of soap residue and lotion, and allowed to dry thoroughly before applying UltraDrape / UltraDrape II to prevent skin irritation and ensure good adhesion.
- Clipping of hair from application site may improve dressing adhesion; shaving is not recommended.
- UltraDrape / UltraDrape II may remain on skin for up to seven days. Replace dressing if it becomes damp, loosened, or visibly soiled.
- For additional information, refer to UltraDrape / UltraDrape II Application Instructions.

UltraDrape Removal Technique

- 1) Gently loosen top edge of dressing.
- 2) Place fingers on top of dressing to stabilize skin or catheter.
- 3) Remove dressing "low and slow" in direction of hair growth, keeping it close to the skin surface and pulled back over itself while stretching the dressing parallel to skin.

IMPORTANT: Removing dressing at an angle may pull at skin, increasing risk of skin injury.

Storage Conditions

Store between 10°C and 30°C, away from direct sunlight.

Disposal

Follow all federal, state and local regulations for non-hazardous waste disposal.

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