


UltraDrape®

REF 34-15, 34-10

EN Instructions for use	JA 取扱説明書
ES Instrucciones de uso	KO 사용법
DE Gebrauchsanweisung	ZH 使用说明
FR Mode d'emploi	AR إرشادات الاستعمال
IT Istruzioni per l'uso	BG Описание
PT Instruções de utilização	HR Opis
DA Brugsanvisning	CS Popis
FI Käyttöohjeet	LV Apraksts
NL Gebruiksaanwijzing	LT Apibūdinimas
SV Bruksanvisning	RO Descriere
EL Οδηγίες Χρήσης	DA opis
HU Használati útmutatás	SV Beskrivning
PO Instrukcja stosowania	TR Açıklama
HE הוראות שימוש	



PARKER LABORATORIES, INC.



Parker Laboratories, Inc.
286 Eldridge Road
Fairfield, NJ 07004 USA
(973) 276-9500
parkerlabs.com



Patent parkerlabs.com/ultradrape



Parker Laboratories BV
Bedrijvenpark Twente 305
7602 KL Almelo
The Netherlands

EN

UltraDrape® and UltraDrape® II UGPIV Barrier and Securement

Instructions for Use:

Description

UltraDrape / UltraDrape II is a uniquely designed, sterile, dual-action barrier and securement dressing for use during Ultrasound-Guided Peripheral Intravenous (UGPIV) procedures. UltraDrape is ASAP (The Association for Safe Aseptic Practice) approved for facilitating Standard - ANTT® (Aseptic Non-Touch Technique). UltraDrape is non-irritating, non-sensitizing and is not made with natural rubber latex. UltraDrape is 5.8 in (14.72 cm) x 3.25 in (8.25 cm). UltraDrape II is 4.68 in (11.9 cm) x 2.63 in (6.7 cm).

Ultrasound gel is applied to removable film layer to prevent contamination of injection site. Convenient adhesive strips on the release film provide added securement of tubing and catheter hub surround. The bifurcated, "stand-alone" design prevents gel from reaching the IV site while enabling seamless catheter insertion and securement.

Indications for Use

UltraDrape / UltraDrape II is a sterile, single-use disposable device, serves as a sterile, viral barrier when used in conjunction with ultrasound transducers to shield patients from the potential transfer of contaminants and biological risk agents that may be introduced during ultrasound scanning and needle-guided procedures. It also serves as a transparent dressing to protect and secure devices to the patient's skin when required. UltraDrape / UltraDrape II is indicated for adults and pediatric patients. For use by healthcare professionals appropriately trained in ultrasound within healthcare settings utilizing ultrasound.

Any serious incident that has occurred in relation to device should be reported to the manufacturer and competent authority of the EU Member State in which user and/or patient is established.



Precautions

- Discard product if package is damaged or evidence of contamination is present.
- Do not re-use.
- Do not re-sterilize.

UltraDrape / UltraDrape II Application Considerations

- Assure skin is clean, free of soap residue and lotion, and allowed to dry thoroughly before applying UltraDrape / UltraDrape II to prevent skin irritation and ensure good adhesion.
- Clipping of hair from application site may improve dressing adhesion; shaving is not recommended.
- UltraDrape / UltraDrape II may remain on skin for up to seven days. Replace dressing if it becomes damp, loosened, or visibly soiled.
- For additional information, refer to UltraDrape / UltraDrape II Application Instructions.

UltraDrape Removal Technique

- 1) Gently loosen top edge of dressing
- 2) Place fingers on top of dressing to stabilize skin or catheter.
- 3) Remove dressing "low and slow" in direction of hair growth, keeping it close to the skin surface and pulled back over itself while stretching the dressing parallel to skin.

IMPORTANT: Removing dressing at an angle may pull at skin, increasing risk of skin injury.

Storage Conditions

Store between 10°C and 30°C, away from direct sunlight.

Disposal

Follow all federal, state and local regulations for non-hazardous waste disposal. Store between 10°C and 30°C, away from direct sunlight.

Disposal

Follow all federal, state, and local regulations for non-hazardous waste disposal.