Aquasonic® 100 Ultrasound Transmission Gel Aquasonic CLEAR® Ultrasound Gel

REF 01-02, 01-08, 01-34, 01-50

REF 03-02, 03-08, 03-34, 03-50, 03-54

- EN** Instructions for use
- ES** Instrucciones de uso
- DE** Gebrauchsanweisung
- FR** Mode d'emploi
- IT** Istruzioni per l'uso
- PT** Instruções de utilização
- DA** Brugsanvisning
- FI** Käyttöohjeet
- NL** Gebruiksaanwijzing
- SV** Bruksanvisning
- EL** Οδηγίες Χρήσης
- HU** Használati útmutatás
- PO** Instrukcja stosowania
- HE** הוראות שימוש
- JA** 取扱口明書
- KO** 사용법
- ZH** 使用说明
- AR** إرشادات الاستعمال
- CZ** Návod k použití
- NO** Bruksanvisning

PARKER LABORATORIES, INC.

Parker Laboratories, Inc.
286 Eldridge Road
Fairfield, NJ 07004 USA
(973) 276-9500
parkerlabs.com

Parker Laboratories BV
Bedrijvenpark Twente 305
7602 KL Almelo
The Netherlands

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

EN

Aquasonic® 100 Ultrasound Transmission Gel
Aquasonic® Clear Ultrasound Gel

Description

Aquasonic 100 Ultrasound Transmission Gel and Aquasonic Clear Ultrasound Gel (Parker Laboratories Ultrasound Gel) are non-sterile, water-soluble and non-staining gels intended for external use only and are used with all commercially available ultrasound devices.

Instructions for Use

In accordance with the facility's hygiene protocol, apply ultrasound gel to the transducer surface and/or the patient's skin, as appropriate.

Avoid direct contact between gel container dispensing tip and any persons or instrumentation, including the ultrasound transducer.

To remove gel, use disposable paper towel, gauze or lint-free cloth.

Discard empty dispenser after use.

Indications for Use

Parker Laboratories Ultrasound Gel can be used on intact, unbroken skin and on all patients in facilities where cross-contamination is of minimal concern. Refill containers and refillable dispensers (product numbers 01-50, 03-50 and 03-54) are for use on low-risk patients and procedures in non-acute care facilities.

Contraindications

Parker Laboratories Ultrasound Gel, in multi-use packaging, is for external use only and should not be used in:

- Any invasive procedures in which a device is passed through the tissue, e.g. needle aspiration, needle localization, tissue biopsy, or other surgical procedures
- Any procedure on or near broken skin, such as fresh surgical site or an open wound
- Patients with severe immunodeficiency or on immunosuppressive therapy
- Neonates or critically ill pediatric patients
- Mucous membranes
- Critically ill patients or patients or those under contact, airborne or droplet transmission-based precautions

Storage Conditions

Store in a cool dry place. Store away from direct sunlight.

Warnings and Precautions

Do not use on patients with a known sensitivity to propylene glycol, as this could cause an allergic reaction. Discard product if package is damaged or evidence of contamination is present.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

Directions for Refilling Dispenser Bottles

Note: Consult your facility policy before refilling from bulk container.

Refillable containers should be washed in hot soapy water, rinsed and thoroughly dried before being refilled. Prior to filling, wipe bottleneck and dispenser tip with isopropyl alcohol or a hospital-grade disinfectant, such as Protex® ULTRA Disinfectant Wipes.

Refill product into corresponding labeled dispenser, and mark with lot number.

Refillable bottles should be marked and discarded according to facility policy; it is recommended bottles be discarded within one month.

Refillable bottles should not be topped off. Prior to refilling, contents should be consumed, or remaining gel should be disposed.

Cracked or damaged bottles should be discarded immediately.

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