# Spectra® 360 Electrode Gel

REF 12-02, 12-08

- EN Instructions for use
- ES Instrucciones de uso
- DE Gebrauchsanweisung
- FR Mode d'emploi
- IT Istruzioni per l'uso
- PT Instruções de utilização
- DA Brugsanvisning
- FI Käyttöohjeet
- NL Gebruiksaanwijzing
- SV Bruksanvisning
- EL Οδηγίες Χρήσης
- HU Használati útmutatás
- PO Instrukcja stosowania
- HE הוראות שימוש
- JA 取扱説明書
- KO 사용법
- ZH 使用说明
- AR إرشادات الاستعمال
- CZ Návod k použití
- NO Bruksanvisning

PARKER LABORATORIES, INC.

Parker Laboratories, Inc.  
286 Eldridge Road  
Fairfield, NJ 07004 USA  
(973) 276-9500  
parkerlabs.com

Parker Laboratories BV  
Bedrijvenpark Twente 305  
7602 KL Almelo  
The Netherlands

MDSS CH GmbH  
Laurenzenvorstadt 61  
5000 Aarau  
Switzerland

EN

Spectra® 360 Electrode Gel

## Instructions for Use:

### Description

Spectra 360 Electrode Gel is a non-sterile, salt-free, chloride free electrically conductive gel. It is non-irritating and bacteriostatic. Spectra 360's non-gritty formula allows for prolonged use without reapplication. Spectra 360's salt-free characteristics make it particularly suitable for long-term applications. Spectra 360 works by wetting the skin, thereby reducing skin resistance.

### Indications for Use

Spectra 360 is recommended for all electro-medical procedures (except defibrillation), including ECG, TENS, monitoring, and can be used on pediatric patients.

### Contraindications

- Spectra 360 is not indicated for defibrillation procedures.

### Warnings

- Do not use on patients with a known history of sensitivity to parabens.

*Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.*

### Precautions

- For external use only on intact, unbroken skin.
- Discard product if package is damaged or evidence of contamination is present.
- To avoid contamination of product, the tube container tips should not come into contact with patient, staff or instruments.

### Directions for Use

Refer to the Instructions For Use provided by the electrode manufacturer of the electrodes, TENS unit, or other device being used for information on how Spectra 360 should be used with their product. Spectra 360 is water soluble, and can be removed using a moistened soft cloth.

### Handling

No special handling required. Spectra 360 is non-flammable and non-corrosive.

### Storage

Store at ambient temperature. Protect product from direct sunlight.

### Disposal

Follow all federal, state, and local regulations for non-hazardous waste disposal.

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